Hammacher
Instrumente

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97 CATALOGUE

ORTHODONTICS

Solingen

Hammacher

Instrumente

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Cleaning, disinfection, drying

Hammacher pliers are suitable for all standard cleaning and disinfection processes like immersion bath, ultrasonic cleaning and thermal disinfection.

The following procedures should always be observed:

When using chemical disinfectants please make sure that they are put together in accordance with the manufacturer's instructions and that the correct concentrations will be used. Please prepare disinfectants freshly every day to avoid evaporation and contamination. It should be ensured that the concentrate has completely dissolved before placing the instruments in the bath. The maximum treatment time and temperature (not above 25 °C / 77 °F) should not be exceeded. Instruments should never be left in the bath overnight or over the weekend !!

Never use any hard objects like metal brushes when cleaning manually, as this could destroy the surface of the instruments and will increase the danger of oxidation. Do not exceed the maximum treatment period when using ultrasonic cleaning equipment.

A suitable non-corrosive solution should be used. Acidic chemicals – especially cement removers – may cause irreparable damages to your instruments. Pliers and instruments with Tungsten Carbide inserts must be disinfected with especially mild chemicals containing corrosion-protection agents !!

The pliers and instruments must be carefully cleaned immediately after the chemical disinfection and cleaning, to remove all chemical residues. They should be rinsed thoroughly in distilled water containing no salt whatsoever. The instruments should be left with joints in the fully opened position !!

We recommend the use of an enzymatic detergent that is pH-neutral and low foaming in case of machine cleaning the instruments. It should be used demineralized water for the final rinsing in the instrument washer.

We recommend the use of a detergent that is neutral or has a low alkaline value between pH7 and 10 in the case of manual cleaning of the instruments

Only high quality distilled water according to EN 285 should be used !! Other types of water may have a too high silicate content which could lead to coloured stains on the instruments during the sterilization process. If rinsing is neglected the following sterilization operations may cause discolouring or rusting of the instruments. Then the instruments must be dried under consideration that the inside of the joints must be dried with special care – we recommend blowing through with compressed air.

Thermo disinfection

Always add a sufficient quantity of the correct chemical when using a thermo disinfectant to neutralize the corrosive effect of the water at high temperatures. Immediately remove instruments from the thermo disinfectant after cleaning and never leave them in damp environment because of the danger of oxidation.

Lubrication

All critical areas on pliers like joints, locks or gliding surfaces for springs should be dried and treated with a protective oil before sterilization. This oil should be permeable by steam, possible to sterilize and temperature-resistant. The treatment with the oil of the Tungsten Carbide inserts reduces the danger of oxidation. We recommend the Hammacher protective oil HSZ 997-12. When using high heat sterilization only lubricate the inner part of the joints and no outside parts to avoid the danger of staining.

Sterilization

Only carefully cleaned and dried pliers and instruments – without any trace of rust – should be placed in the sterilizer or autoclave bearing in mind that instruments should always be left with joints in the opened position. If rusty pliers and instruments are sterilized with clean uncorroded ones, a rust film may appear which can seriously damage previously intact material. Special care is required when sterilizing chromium-plated pliers and instruments – what can be made but it is not recommended because the joints have no chrome-plating and therefore can corrode very quickly. For this reason Hammacher pliers and instruments for the treatment area are made of stainless steel.

Pliers and instruments should not be stored for long periods of time in a damp environment. Therefore the door of the sterilizer should be opened immediately after the sterilization process. If necessary the pliers and instruments can be dried with sterile cloths. When foil bags are being used, make sure that they are correctly positioned in the sterilizer and that they contain no residual moisture.

Sterilization at temperatures exceeding 180 °C / 356 °F may harm the hardness of the tips and cause discoloration of the pliers and instruments. The cutting edges may also become dull more quickly and there may be an increased danger of oxidation

Use of new pliers and instruments

Made of carefully hardened stainless steel and processed by experienced and skilled workmen on modern machinery our pliers achieve the highest quality.

Before using new pliers and instruments they must be thoroughly cleaned and subjected to the entire process of preparation. Compliance with the above instructions will preserve the effectiveness and prolong the lifetime of your Hammacher pliers and instruments.

The guarantee does not apply if the equipment has been deliberately modified or used/cleaned in an unprofessional manner. Improper use or incorrect sterilization including cleaning, processing and lubrication as well as negligence and wear and tear void the guarantee.

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herewith confirms that the company

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has introduced, applies and maintains a Quality Management System in the area of:

**Purchase, manufacture, final inspection, storage and
distribution of reusable dental and surgical instruments**


Distribution of instruments and disposables

The compliance of the Quality Management System with the requirements of the below mentioned standard was verified by an audit:

EN ISO 13485

The license of certification is subject to surveillance by MEDCERT.

Report No.: 3141FS14F
Process No.: QS – 3141
Certificate No.: 3141GB438181030

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

Instruments made of stainless steel are rust proof but not rust-free like the Wironit instruments.

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Titanium instruments are 40% lighter than instruments made of stainless steel - they are absolutely rust-proof and suitable for all common types of cleaning and sterilization.

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Wironit instruments are made of a special alloy and show the following advantages:

- ✓ rust resistant
- ✓ anti-magnetic
- ✓ acid resistant

They are beautifully designed and of unlimited durability.

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All products marked with CE are medical devices.
All products not marked with CE are no medical devices.

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Activating lever	35	Corrosion and condense eraser	48
Adams mod. Sling, clasp and wire bending pliers	61, 62	Crown and collar scissors	26
Adams mod. Wire bending pliers	59	Cutter, Distal end	37, 75, 76
Adams Sling, clasp and wire bending pliers	61	Cutter, Hard wire	78
Aderer Wire and clasp bending pliers	57, 74	Cutter, Ligature	38, 39, 77, 78
Angle Arch forming pliers	68	De La Rosa Arch contouring pliers	65
Angle mod. Sling, clasp and wire bending pliers	61	Diagnostic set	22
Angle/Bird Beak Sling, clasp and wire bending pliers	61, 62	Digital pocket caliper	13
Angle/Tweed Arch forming pliers	68	Distal end vender	29
Arch contouring pliers	65	Distal end cutter	37, 75, 76
Arch forming pliers	68	Double action side nipper	72
Arch tucker	35	Easy Lift	43
ARQUATUS Diagnostic set	22	Eraser, Corrosion and condense	48
ARQUATUS Mouth mirror handle	22	Examination mirror	20
ARQUATUS Probe	22	Flat nose pliers	69, 70
ARQUATUS Tweezers College	22	Foil and splints scissors	26, 42
Artery forceps	34	Forceps, Artery	34
Band / bracket positioner	29	Halsted-Mosquito Artery forceps	34
Band pinshing and ligature locking pliers	54	Hammacher/Aderer Wire and clasp bending pliers	36, 58
Band pusher	29	Hammacher/Oliver Band removing pliers	54
Band removing pliers	54, 55	Handle, Mouth mirror	20, 22
Bar bending pliers, Wire and	60, 74	Hard wire cutter	78
Beebee Crown and collar scissors	26	Heidemann Spatula	30
Beerendonk Sliding caliper	15	Holder, Needle	32, 33
Bending instrument, Wire	29	How Band pinshing and ligature locking pliers	54
Bending pliers, Loop	63	Impression trays	2 - 10
Bending pliers, Pin	73	Inserting combination instrument	30
Bending pliers, Sling, clasp and wire	61, 62	Instrument for removing adhesive residue	29
Bending pliers, Wire	59	Instrument, Inserting combination	30
Bending pliers, Wire and bar	60, 74	Instrument, Ligature	28, 33
Bending pliers, Wire and clasp	36, 57, 58	Instrument, Wire bending	29
Bite stick	29	Instruments holder	44 - 47
Box, Orthodontic	48	Interstice width caliper	11
Box, Retainer	43	Intraoral mirror	18
Box, Wire dispenser	50	Jarabak Wire and bar bending pliers	60
Bracket / band positioner	29	Jarinda Wire and bar bending pliers	60
Bracket placing tweezers	23, 24	Kalthoff Pliers	71
Bracket placing tweezers Schrems	25	Korkhaus Case for kit	12
Bracket positioner	31	Korkhaus cm-disc	12
Bracket removing pliers	35, 56	Korkhaus mm-disc	12
Caliper, Digital pocket	13	Korkhaus Orthometer	12
Caliper, Interstice width	11	Korkhaus Orthometer Kit	12
Caliper, Sliding	15	Korkhaus Space divider with socket	12
Case for kit Korkhaus	12	Ligature cutter	38, 39, 77, 78
Cheek and lip retractor	17	Ligature instrument	28, 33
Cheek retractor	16	Ligature locking pliers, Band pinshing and	54
Cheek retractor and mouth mirror, Combined	19	Ligature pliers, Wire seizing and	39, 51, 52
Clamp	34	Ligature tucker	28
Clasp and wire bending pliers, Sling,	61, 62	Lingual servant	35
Clasp bending pliers, Wire and	36, 57, 58	Lip and cheek retractor	17
Closing loop pliers	63	London-College Tweezers	25
cm-disc Korkhaus	12	Loop bending pliers	63
College Tweezers	22, 25	Loop pliers, Closing	63
Combination instrument, Inserting	30	Malotky Wire and banr bending pliers	74
Combined mouth mirror and cheek retractor	19	Mathieu-Kocher Needle holder	32
Condense and corrosion eraser	48	Measuring probe	22
Contouring pliers, Arch	65	Measuring ruler	11
Contrastor, Palatal	17	Micro-Halsted Artery forceps	34

Mirror, Examination	20	Punch Slot Machine	41
Mirror, Intraoral	18	Punch The GUN I	40
Mirror, Mouth	20, 21	Pusher, Band	29
Mirror, Parallelometer	20	Refiner	42
Mirror, Photo, from glass	19	Removing adhesive residue instrument	29
Mirror, Photo, from metal	18	Removing pliers, Bracket	35, 56
mm-dics Korkhaus	12	Retainer box	43
Mouth mirror	20, 21	Retractor, Cheek and lip	17
Mouth mirror and cheek retractor, Combined	19	Retractor, Cheek	16
Mouth mirror handle	20, 22	Ruler, Measuring	11
Münchner Modell Sliding caliper	15	Schrems Bracket placing tweezers	25
Nance Closing loop pliers	63	Scissors for foils and splints	26, 42
Needle holder	32, 33	Scissors, Crown and collar	26, 27
Nice End Pliers	53	Separating pliers	34
Oil pen	48	Set, Orthodontic starter	79 - 83
Oliver Band removing pliers	54	Side nipper	72
Optical pliers	70	Side nipper, Double action	72
Orthodontic box	48	Sliding caliper	15
Orthodontic probe	30	Slot Machine Punch	41
Orthodontic starter package	79 - 83	Space divider with socket Korkhaus	12
Orthometer Kit Korkhaus	12	Stand for pliers	49, 50
Orthometer Korkhaus	12	Stand for tweezers	50
Palatal contrastor	17	Starter package, Orthodontic	79 - 83
Parallelometer mirror	20	Step pliers	64
Pen, Oil	48	Sterilization cassette	44 - 47
Photo mirror from glass	19	Stop pliers V-bend	64
Photo mirror from metal	18	The GUN I Punch	40
Pin bending pliers	73	Tip protection	48
Placing tweezers, Bracket	23, 24	Torquing pliers	66, 67
Pliers Nice End	53	Trays, Impression	2 - 10
Pliers stand	49, 50	Tucker, Ligature	28
Pliers V-bend, Stop	64	Tweezers College	22, 25
Pliers, Arch contouring	65	Tweezers London-College	25
Pliers, Arch forming	68	Tweezers stand	50
Pliers, Band pinshing and ligature locking	54	Tweezers, Bracket placing	23, 24
Pliers, Band removing	54, 55	Twister	29
Pliers, Bracket removing	35, 56	Universal Crown and collar scissors	27
Pliers, Closing loop	63	Universal pliers	73
Pliers, Flat nose	69, 70	Vestibulum Cheek and lip retractor	17
Pliers, Loop bending	63	Waldsachs Universal pliers	73
Pliers, Optical	70	Weingart Wire seizing and ligature pliers	39, 51, 52
Pliers, Pin bending	73	Wilhelmy Activating lever	35
Pliers, Separating	34	Wilhelmy Arch tucker	35
Pliers, Sling, clasp and wire bending	61, 62	Wilhelmy Bracket removing pliers	35
Pliers, Step	64	Wilhelmy Lingual servant	35
Pliers, Torquing	66, 67	Wilhelmy/Hammacher Clamp	34
Pliers, Universal	73	Wire and bar bending pliers	60, 74
Pliers, Wire and bar bending	60, 74	Wire bending instrument	29
Pliers, Wire and clasp bending	36, 57, 58	Wire bending pliers	59
Pliers, Wire bending	59	Wire bending pliers, Sling, clasp and	61, 62
Pliers, Wire seizing and ligature	39, 51, 52	Wire dispenser box	50
Pocket caliper, Digital	13	Wire seizing and ligature pliers	39, 51, 52
Positioner, Band / bracket	29	Young Loop bending pliers	63
Positioner, Bracket	31	Zürcher Modell Sliding caliper	15
Probe	22, 30		
Probe, Measuring	22		
Probe, Orthodontic	30		
Prof. Dr. Günay Cheek retractor	16		

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
HWL 374-18
(OK-18 XS)

CE

L =	groß / Kind large / child
M =	mittel / Kind medium / child
S =	klein / Kind small / child
XS =	extra-klein / Kind extra small / child

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Set with 8 children sized impression trays, consisting of:
4 upper jaws: **HWL 371-15 to HWL 374-18**
4 lower jaws: **HWL 376-15 to HWL 379-18**

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Set with 8 children sized impression trays, perforated, consisting of:
4 upper jaws: **HWL 381-15 to HWL 384-18**
4 lower jaws: **HWL 386-15 to HWL 389-18**

Regular

**Oberkiefer
Upper jaws**

HWL 301-01
(OK-1 XL BS5 US-U11)

CE

HWL 302-02
(OK-2 L BS3 US-U12)

CE

HWL 303-03
(OK-3 M BS2A US-U13)

CE

L =	extra groß extra large
L =	groß large
M =	mittel medium
S =	klein small
XS =	extra klein extra small

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Regular

**Unterkiefer
Lower jaws**

XL =	extra groß extra large
L =	groß large
M =	mittel medium
S =	klein small
XS =	extra klein extra small

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Set with 10 impression trays, consisting of:
5 upper jaws: **HWL 301-01 to HWL 305-05**
5 lower jaws: **HWL 311-01 to HWL 315-05**

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"Wiro -lock" Impression Trays

- ✓ "Wiro -lock" Impression Trays are made of the high value, corrosion proof WIRONIT Steel.

An absolute high-quality-product !

- ✓ These "Wiro -lock" Impression Trays dispose of a retention edge lock. So they are able to hold the i mpression material without using any adhesives.
- ✓ "Wiro -lock" Impression Trays might be cleaned inside a supersonic bath and of course they are sterilizable by any usual methods.

perforiert
perforated

Oberkiefer
Upper jaws

HWL 306-01
(OK-1P XL PS6 US-7)

CE

HWL 307-02
(OK-2P L PS5 US-5)

CE

HWL 308-03
(OK-3P M PS4 US-4)

CE

HWL 309-04
(OK-4P S PS3 US-1)

CE

XL =	extra groß extra large
L =	groß large
M =	mittel medium
S =	klein small

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Set with 8 impression trays, consisting of:

4 upper jaws: **HWL 306-01 to HWL 309-04**
4 lower jaws: **HWL 316-01 to HWL 319-04**

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
Interstice width caliper

for easy and precise measurement under usage of extension plates and activators.

Measuring range 4mm to 80mm.
for comparative measuring

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Measuring ruler

measurement area of 80 mm

(with mm-scale on the upper and lower border)
The dental ruler for the registration, examination and documentation.
for comparative measuring

HSL 248-70
110 mm / 4-5/16 inch
sterilisierbar / sterilizable

Messlineal aus rostfreiem Stahl

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Measuring ruler from stainless steel

Measuring range: 0 - 60 mm

for comparative measuring

HSL 252-15
15 cm / 6 inch
sterilisierbar / sterilizable

HSL 252-20
20 cm / 8 inch

Messlineal aus rostfreiem Metall

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Measuring ruler of from stainless steel

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
Orthometer Kit

Based on the Pont-index, modified according to Korkhaus, this kit determines the normal numbers as an indication for building a regular dental arch.

Consisting of:

Orthometer
Space divider
cm- and mm-disc
Case

All parts are available separately.

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
Space divider with socket

Brass chromated, with steel tip
only disinfectable

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
- ✓ with narrow shaped measuring jaws
- ✓ with cross peaks and a rod to measure the depth
- ✓ from stainless steel, hardened and tempered
- ✓ jumbo display, reading 0,01mm / .0005"
- ✓ switch on / off, DIN 862, metal casing, incl. wallet, without roller
- ✓ measuring range: 100 mm, jaws: 30 mm

SERVICE

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We consider **SERVICE** one of the main aspects of our philosophy!

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


extraoral

HSL 247-50
Beerendonk
125 mm / 4 3/4 inch
sterilisierbar / sterilizable

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Sliding caliper
Measuring range 0-80mm

HSL 247-51
Münchner Modell
125 mm / 4 3/4 inch
sterilisierbar / sterilizable

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Sliding caliper
Measuring range 0-80mm

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Sliding caliper
Measuring range 0-80mm.

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


Set of 2 palatal contrastors (1 wide + 1 narrow)

- From black anodised aluminium.
- The contrastors are suitable for normal and digital photography.
- They give contrast to your images and eliminate those elements not required in the picture.

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Lip and cheek retracor, pair

for adults, made of autoclavable polycarbonat. Easy to snap into place. Allows a clear view and working area. Comfortable for the patient. Non-reflective, ideal for photography. For sterilization we recommend the safety programme of 121°C - please avoid metal instrument contact.

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
Cheek and lip retractor / Vestibule-hook

for the entire vestibule of the upper and lower jaws

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
Lip and cheek retracor, pair

for children, made of autoclavable polycarbonat. Easy to snap into place. Allows a clear view and working area. Comfortable for the patient. Non-reflective, ideal for photography. For sterilization we recommend the safety programme of 121°C - please avoid metal instrument contact.

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Photo mirror from metal with elaborate rhodium coating
ideal for distortion-free intraoral images
durable and robust
sterilizable at 134° C
with removable silicone handle

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Double-sided photo mirror from glass, 3mm thick

- ✓ The mirrors can be treated in a hot air sterilizer up to 180° C, in an autoclave up to 130° C or in a thermo disinfectant up to 105° C.
- ✓ The removable silicone handle is autoclavable.
- ✓ The coating is made from an extremely glossy float glass.
- ✓ The mirrors provide almost 100% reflection colour fidelity and are corrosion resistant.

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ULTRA FS – the brightest and sharpest mirror in the world – the new unique design of the reflective layers offers an unsurpassed bright reflection. 1 pack = 6 pcs.

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The ULTRA FS quality features are:

- ✓ Highest reflection for a brilliant brightness
- ✓ Superior colour accuracy without any colour distortion
- ✓ Exactly sharp mirror image
- ✓ Scratch-resistant
- ✓ Optimized holding of the cheek

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


Mouth mirror measuring-handle
hollow, hexagonal handle with etched
80mm-scale, one end is oblique.
For measuring the levels of the
orthodontic indication groups.

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„HAMSO“ Front Mouth Mirrors

- Front surface, plane (1 package = 12 pieces)
- Very bright and clear picture
- Reflection over 90 %
- Colour fidelity
- Highly scratch-resistant
- Suitable for all usual kinds of disinfection and sterilisation
- Temperature-resistant up to 200°C/390°F

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
Mirrors from glass rhodium coated

- ✓ PP plastic handle reinforced with glass fiber

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RELAX | Mouth Mirrors | 1 pack = 10 pcs. | different colors (arbitrarily) CE

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RELAX | Mouth Mirrors

The benefits of a front surface mirror in combination with the "lite" construction of the fibre optic provide you and your patient with the highest level of comfort during the treatment. The instrument is autoclavable, acid-resistant and suitable for the thermal disinfectant.

- ✓ **Pleasantly light** - for a relaxed work
- ✓ **Optimized holding** - thanks to the wide handle with hollows
- ✓ **Quality resistant**, colour stable, acid resistant and resistant to plaque indicators
- ✓ **Patient friendly** - thanks to the rounded shapes
- ✓ **Hygienic** - thanks to the smooth transition from mirror to frame
- ✓ **silk mat**, non-reflecting surface

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Probes

Measuring probes, single
ended, with 3 mm
graduations, for the
measurement of the
orthodontic indication groups

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Placing tweezers
for an easy placing of the
brackets – especially for bonding
tubes

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**Double angled placing tweezers, crossed, with
cleaning scaler**

This placing tweezers is particularly suitable for bonding of
the brackets in the molar range. With its short and narrow
working ends you can seize under the wings of the
brackets and place the brackets exactly without a slipping
of the bracket when releasing. Because of the double angle
the sight is less limited.
The rear delicate working end is particularly suitable for
opening and closing of selfligaturing brackets, for removing
the adhesive residues and for placing the brackets.

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**Tweezers, crossed,
self retaining / spring action**
Bonding tweezers with delicate working
ends for an easy placing of the brackets
and especially of the molar bondable
tubes.

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
Double angled placing tweezers, crossed,
This placing tweezers is particularly suitable for
bonding of the brackets in the molar range. With
its short and narrow working ends you can seize
under the wings of the brackets and place the
brackets exactly without a slipping of the bracket
when releasing. Because of the double angle the
sight is less limited.

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
Bonded bracket placing tweezers
spring action. The joint serves for the angulation of edgewise-brackets.

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Bracket placing tweezers, small
suitable for lingual treatment.

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
Indirect bonding tweezers
for lingual techniques

For bonding lingual brackets which were placed by transfer flaps. Therefore an O-ring as "handle" is required. Because of the bevelled working end of the tweezers the area of view is not limited and the exact placing of the transfer flap is guaranteed.

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Tweezers for placing molar brackets acc. to Dr. Dr. Schrems

With this precision tweezers you can place brackets to all teeth of the upper and lower jaws (from the front teeth of the upper jaws to the molars – including the 2nd molars) with an optimum of visibility. After opening the tweezers there are no unintended changes in the positioning brackets. After placing the brackets with these tweezers there is only the fine tuning with the usual instruments necessary.

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
Scissors for foils and splints, curved
with micro saw edge, **stainless**.
Suitable for all types of foils – for Aligner splints.

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Scissors for foils and splints, curved
with micro saw edge, chromium plated.
Only disinfectable!
Suitable for all types of foils – for Aligner splints.

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Ligature tucker
with band pusher
CE

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
Band- / bracket positioning gauge

Easy seating of the bracket on the teeth with exactly the right distance between bracket slot and incisal edge. Four different heights for bracket positioning: 3,5mm, 4mm, 4,5mm, 5mm.

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Band seater
Autoclavable bite stick for adjusting and inserting bands

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Wire bending instrument

This instrument facilitates the arch wire bending behind the tubes. The wire can be bended in every direction without damaging the wire.

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Inserting combination instrument

With the fork on one side of the instrument you can seize the arch distal and mesial and you can easily put it into the slot of the brackets – especially with rotated teeth. At the other end it is a traditional ligature carrier to hide the wire ligatures.

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Orthodontic probes with a slightly thicker handle
for the relief of the wrist. Very short probes for easy removal of ligatures and elastics.

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Orthodontic probes with PPSU handle
for the relief of the wrist. Very short probes for easy removal of ligatures and elastics.

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Probe and explorer, Wironit
single ended with octagonal handle

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Delicate needle holders
for ligating dentalastics

CE

HSD 009-14
Mathieu-Kocher
140 mm / 5-1/2 inch
sterilisierbar / sterilizable

Nadelhalter
Needle holder

CE

HSD 170-14
Tungsten Carbide Inserts (TC)
Mathieu-Kocher
140 mm / 5-1/2 inch
sterilisierbar / sterilizable

Nadelhalter, schlankes Maul
Needleholder, with slender jaws

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Tying and Elastomeric plier

Inserted Mathieu style.
Very thin tips for a clear view. All stainless steel. Very thin tips allow you to grasp ligatures or elastomerics with maximum visibility. Spring action in the center of the handle allows smooth travel without any sticking or dragging during ligation.
For wire ligatures from .008" bis .014"
The ideal instrument for ligation!

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
Tying and Elastomeric pliers, with very fine tips

Inserted Mathieu style.
Very thin tips for a clear view. All stainless steel. Very thin tips allow you to grasp ligatures or elastomerics with maximum visibility. Spring action in the center of the handle allow smooth movement without any sticking or dragging during ligation.

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
Separating pliers
for expanding Dentalastics before positioning.

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
Bracket removing pliers

These bracket removing pliers are especially destined for the lingual orthodontia. The angled and very delicate working ends make the removing of the brackets much easier. The working ends will be put directly at the base of the lingual bracket which allows a gentle removal. The long handles of the pliers also ease the removal of the brackets 6 and 7.

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The arch tucker serves for setting and bending of arches.

CE

HSL 803-00
170 mm / 6 11/16 inch
Wilhelmy
Aktivierhebel .016" x .022"
sterilisierbar / sterilizable

Der Aktivierhebel (Finisher) dient zum intraoralen Aktivieren 16x22 SS (Artistik).

The activating lever (Finisher) serves for intra-oral activating 16x22 SS (Artistik).

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The Lingual Servant serves for ligating in and out (Alastics).

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Distal end cutter with safety hold

- cuts the wire directly behind the tubes or brackets – the cut wire is being held by the safety hold
- min. Ø □ 0,36 mm / .014"
- max. Ø □ 0,43 mm / .017"

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Distal end cutter with safety hold

- cuts the wire directly behind the tubes or brackets – the cut wire is being held by the safety hold
- min. Ø □ 0,43 mm / .017"
- max. Ø □ 0,48 mm / .019"

HSL 2827-16
Tungsten Carbide Inserts (TC)
165 mm / 6-1/2 inch
sterilisierbar / sterilizable

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Ligature cutter

with **outside** cutting edge

With the pointed working end of the angled ligature cutter it is possible to separate the wire directly at the bracket. The angle makes the work much easier as the cheeks and labial angles have not to be pulled backwards unnecessarily.

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
This punch is indispensable for the invisible splint therapy. It is used for the easy and reproducible removal of superextended splint borders without the production of grinding spots and sharp edges. The result is a clearly reduced expenditure of time compared to the use of rotary instruments.

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
This punch is suitable for the use of inter- and intra-maxillary elastic strips in connection with deep-drawing splints. A complex positioning of buttons is no more necessary. The special cut-out prevents a slipping out of the elastic strips and a cracking of the splints.

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Hammacher **Refiner** according to Dr. Drechsler
The aimed use of the **Refiner** allows to support individual tooth movements. By creating a point shaped elevation in the aligner a selective activation takes place which increases the pressure on individual teeth. These cold deformation of the aligner in the refinement phase provides great help in the upgrading of the treatment results with the aligner therapy.

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Scissors for foils and splints, curved, with micro saw edge, **stainless**. Suitable for all types of foils – for Aligner splints.

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
Hammacher **EasyLift** according to Dr. Drechsler
EasyLift makes it particularly easy for patients to take out their aligner. Because of special dental morphology or due to appropriate numbers and shapes of glued attachments this is often difficult without help. The delicate rounded hook of the EasyLift can be placed comfortably at the edge of the aligner so that a safe distance is always possible. EasyLift – with its high-quality steel handle – is at the same time very handy so that it can be taken discreetly in the aligner box.

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

Instrument holder

stainless steel and silicone – for 10 instruments – can be used with any known disinfection solution – resistant to any sterilisation and cleaning process

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
Instrument holder

stainless steel and silicone – for 22 instruments – can be used with any known disinfection solution

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


✓ **EASY-TRAY** has a locking arm rather than a loose clip. The arm is hinged and swings up and over the instruments, giving you full and easy access to your instruments. Once used, swing the arm back thus locking it in place, keeping your instruments safe and your staff protected whilst processing for readiness again.

Manufactured from a high quality stainless steel and silicone the **Easy-Tray** is getting by electro-polishing a bright, shiny, smooth and easy to clean surface. **EASY-TRAY** is designed for the use in ultrasonic baths and washer disinfectors with all known solutions and autoclaves to 134°C.

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Advantages:


- ✓ The internal silicone strips put no pressure on the mirrors which allow complete penetration by the washing and disinfection solutions.
- ✓ These cassettes are made of high performance materials (stainless steel and silicone).
- ✓ These cassettes have been designed for use with all disinfection and cleaning solutions such as ultrasonic baths or washer disinfectors.
- ✓ Autoclavable up to 134°C
- ✓ Because they are very thin (1,3cm) these cassettes take up very little room in any process.

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
Multifunctional tray

This multifunctional tray was designed as a normtray insert, consisting of two instrument supports placed one above the other. The normtray insert can be used as a washtrey in the thermal disinfecter and is autoclavable. The plastic version of the tray protects the delicate instrument tips from damage. The instruments lie absolutely free which guarantees an optimal cleaning result. Silicone strips are deliberately omitted in these trays.

For secure fixation and storage of instruments during daily work as well as for gentle mechanical cleaning, disinfection and sterilization. Available in 3 colours and 3 sizes

- ✓ fit in aluminum and stainless steel normtrays
- ✓ suitable for instruments up to 172,5 mm overall length
- ✓ protection against cuts and punctures
- ✓ form stable
- ✓ due to the large perforation / openings a good water-jet permeability is possible (almost without rinsing shadow)
- ✓ high quality plastic especially for medical purposes
- ✓ protection and organisation of the instruments
- ✓ the diagonally arranged meshes facilitate the identification of the parallel placed instruments
- ✓ simple and comfortable handling
- ✓ colour coded for the respective practical application

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
HSM 241-03
rot / red

CE


Tray, small
Outer dimension: 180 x 92 x 34 mm
Inner dimension: 172 x 85 x 29 mm

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Oil pen

- ✓ synthetic maintenance oil
- ✓ vapour permeable and high temperature stability up to 260 ° C
- ✓ odourless, transparent and free of silicone
- ✓ does not leave any brown stains on the instruments
- ✓ precise dosage by the safety / anti-drip valve

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




Wire seizing and ligature pliers
 Holds round, square and rectangular wires
 up to 0,56 mm x 0,70 mm / .022" x .028"

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Wire seizing and ligature pliers
 For holding of wires round, square and rectangular until max. 0,56 mm x 0,71 mm / .022" x .028"

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Universal pliers with small working ends
 • for hard wire up to max. 0,5 mm / .020" Ø

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Hammacher Nice-End-Pliers according to Dr. Banach

Wire arches bending pliers (Cinch-Back-Pliers) suitable for cross sections of wire up to 0,4 mm / .016 "

Nice-End – the pliers to bend highly elastic wire arches with minimal cross section. Up to now impossible – these pliers now bend also thinnest highly elastic arches exactly perpendicular and prevents the typical irritations of the buccal mucosa by stinging arch ends and at the same time secures the wire from slipping through.

4 mm

HSL 223-14
no. 110
How
 140 mm / 5-1/2 inch
 sterilisierbar / sterilizable

Band- und Ligaturenzange
 gerieft

Band pinshing and ligature
locking pliers serrated

CE

4 mm

HSL 226-14
no. 111
How
 140 mm / 5-1/2 inch
 sterilisierbar / sterilizable

Band- und Ligaturenzange
 gerieft

Band pinshing and ligature
locking pliers serrated

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
Wire and clasp bending pliers
• for wire of thicknesses 0,25 mm /
.010" and 0,3 mm / .012" Ø

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Wire and clasp bending pliers
• for hard wire until max.
1,1 mm / .043" Ø

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Set consisting of 2 pliers HSL 4192-13 A and HSL 4192-13 B.

The pliers can also be ordered individually.

HSL 4192-12 A
Hammacher/Aderer
 120 mm / 4-3/4 inch
 sterilisierbar / sterilizable

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Wire bending pliers

- with two grooves in front of the tip – with additional comfort groove for intrusion
- especially for triangular clasps, Adams clasps and springs of the protrusion
- for hard wire until max. 0,9 mm / .035" Ø

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
Wire bending pliers

- with one groove in front of the tip
- especially for triangular clasps, Adams clasps and springs of the protrusion
- for hard wire until max. 0,9 mm / .035"

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Wire and bar bending pliers
• for hard wire until max.
0,5 mm / .020" Ø

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Light wire bending pliers

Fine serrations on square beak and three precision grooves allow accurate bending & closing of loops with a secure grip. Especially good when used on light wire sizes.

Bends round, square and rectangular wires up to 0,5 mm / .020".

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Wire and bar bending pliers

- universally applicable, allows easy and comfortable working in the laboratory
- for all commercially available orthodontic wires until max. 0,8 mm / .031" (hard and spring hard)

HSL 401-13
no. 139
Angle / Bird Beak
130 mm / 5-1/8 inch
sterilisierbar / sterilizable

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Sling, clasps and wire bending pliers
for hard wire up to 0,7 mm / .028" Ø max.

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Wire bending pliers

Excellent for all types of wire bending in the laboratory.
 Stubby beaks - round on square.
 Bends round, rectangular and square wires up to 0,76 mm / .030".

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Bird beak / serrated

The non-slip grooves on the square jaw part offer you maximum control over the wire. With these pliers you do not put pressure on the wire, but hold and fix it with the jaw ends. The entire side of the jaws is serrated and therefore these pliers are ideal for the bending of retaining clasps and wires which have the tendency to slip.
 For bending of round, square and rectangular wire up to max. 0.76 mm / .030".

HSL 406-14
Tungsten Carbide Inserts (TC)
Adams modif.
 125 mm / 4-7/8 inch
 sterilisierbar / sterilizable

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Sling, clasps and wire bending pliers
 for hard wire up to 0,7 mm / .028" Ø max.

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Closing loop pliers
for hard wire until max.
0,6 mm / .024" Ø

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Stop pliers V-bend

- for placing a 1 mm stop
- for round, square and rectangular wire until max. 0,56 mm / .022"

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Step pliers

- automatic bending of 0,5 mm steps
- for round, square and rectangular wire until max. 0,53 x 0,64 mm / .021 x .025"

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Step pliers

- automatic bending of 0,75 mm steps
- for round, square and rectangular wire until max. 0,53 x 0,64 mm / .021 x .025"

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Step pliers

- automatic bending of 1,0 mm steps
- for round, square and rectangular wire until max. 0,53 x 0,64 mm / .021 x .025"

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
De La Rosa arch contouring pliers, non-grooved
Non-grooved pliers allow contouring the archwire at any position on the instrument.
Bends round, square and rectangular wires up to 0,5 mm / .020".

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
Set consisting of 2 pliers HSL 4246-13 A and HSL 4246-13 B.
The pliers can also be ordered individually.

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Single tooth torque set for anterior mandibles for 3,5 mm bracket width

- The orthodontist can set the required torque directly during the removal of the arch without twisting the remaining wire.
- A displacement resp. debonding of the brackets is not necessary.
- Allows a time-saving work.
- Only for square wires

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Torquing Plier Set for 4,5mm bracket width

Torques archwires accurately without distortion to the remaining wire.

No need for different sized torquing keys.

- for square and rectangular wire until max. 0,56 x 0,70 mm / .022 x .028"

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NEW

All torque pliers now with
TORQUE INDICATOR

- ✓ With the left and right sides lasered **TORQUE INDICATOR** you can estimate the torque value of 0° to 15° marked in the wire arches.
- ✓ The **TORQUE INDICATOR** increases the precision in the practice.

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NEW

All torque pliers now with
TORQUE INDICATOR

- ✓ With the left and right sides lasered **TORQUE INDICATOR** you can estimate the torque value of 0° to 15° marked in the wire arches.

- ✓ The **TORQUE INDICATOR** increases the precision in the practice.

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Arch forming pliers
for hard wire up to 0,7 mm / .028"

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Tweed arch pliers / 2 Step pliers

Shorter jaws give better control and grip.
Excellent for detailing an arch wire.

Bends round, square and rectangular wires up
to 0,56 mm x 0,7 mm / .020 x .028".

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stainless, serrated
for wire until 0,7 mm / .028" max.

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Flat nose pliers, pointed, narrow, serrated
for hard wire until 0,5 mm / .020" Ø max.

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Flat nose pliers with guide notch,
serrated
for hard wire until 0,7 mm / .028" max.

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Flat-nose pliers, tapered tips
(**Dolphin pliers**), serrated
for wire until 0,7 mm / .028" max.

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Flat nose pliers, serrated, more delicate
 for hard wire up to 0,7 mm / .028" Ø max.

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Biegen einer Frontalfeder
Bending an anterior spring

Retention der Drahtenden
Retention of wire ends

Draht abmessen und trennen
Measuring and cutting of wire
to length

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





Universal pliers

Bends: hard wire up to 1,0 mm / .039" Ø max.
hard profiled wires up to 1,3 x 0,65 mm / .051" x .026" max.
Cuts: hard wire up to 0,7 mm / .028" Ø max.

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Wire and bar bending pliers
 • for hard wire until max. 1,5 mm / .059" Ø

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Wire and bar bending pliers
 • for hard profiled wire until max. 3,0 x 1,5 mm / .059"

HSL 3315-14
Tungsten Carbide Inserts (TC)
 140 mm / 5-1/2 inch
 sterilisierbar / sterilizable

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- min. Ø □ 0,31 mm / .012"
- max. Ø □ 0,64 mm / .035"

Distal end cutter with safety hold, TC

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Ligature cutter

- the delicate form allows an easy separation of the wire ligatures between the bracket wings
- for wire ligatures from min. 0,25 mm / .010" until max. 0,4 mm / .016" Ø

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Ligature cutter

- the small size head is 15 degrees angled for a better visibility and easy access
- for wire ligatures from min. 0,25 mm / .010" until max. 0,4 mm / .016" Ø

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

Instrument holder with cover
stainless steel and silicone – for 10 instruments – sterilizable

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Instrument holder
stainless steel and silicone – for 22 instruments –
can be used with any known disinfection solution

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Wire seizing and ligature pliers
For holding of wires round, square and rectangular
until max. 0,56 mm x 0,71 mm / .022" x .028"

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Bird beak / serrated

The non-slip grooves on the square jaw part offer you maximum control over the wire. With these pliers you do not put pressure on the wire, but hold and fix it with the jaw ends. The entire side of the jaws is serrated and therefore these pliers are ideal for the bending of retaining clasps and wires which have the tendency to slip.

For bending of round, square and rectangular wire up to max. 0.76 mm / .030".

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Hard wire cutter

- with a 15° angle
- for an easier intra-oral access to the posterior area
- for hard wire from min. 0.6 mm / .024" until max. 0.8 mm / .031" Ø

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Reparatur service

The services of Karl Hammacher GmbH also include a repair service, for example, comprising the following work:

- Replacement of hard-metal inserts
- Sharpening and resharpening of scissors, scalpels and other cutting instruments
- Complete reworking of older instruments such as repolishing, straightening, ensuring ease of movement, derusting
- Surface treatment of instruments such as: dimming, polishing, chrome-plating, gold-plating
- Replacement of screws, pins or rivets

Moreover we produce specially for you according to your pattern, sketch or suchlike.

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© copyright: Chantal Design '92 (NL)
In Auftrag von: _ Karl Hammacher GmbH, Solingen